

DISTRIBUTED BY

INSTRUMED SURGICAL

2180 Dunwin Drive, Units 5 & 6
Mississauga, ON L5L 5M8
905-820-0902 (phone)
800-667-5653 (toll free)
orderdesk@instrumedurgical.com
cust.service@instrumedurgical.com




EverLift[™]
Submucosal Lifting Agent

THE LIFT YOU'VE BEEN WAITING FOR

EverLift™

Submucosal Lifting Agent



INDICATION FOR USE


EverLift Submucosal Lifting Agent is indicated for use in gastrointestinal endoscopic procedures for submucosal lift of polyps, adenomas, early-stage cancers, or other gastrointestinal lesions prior to excision with a snare or other appropriate endoscopic device.

SUPERIOR CUSHION


CELLULOSE BASED & HYPERTONIC

EverLift is cellulose based enabling the solution to stay in place longer, as it is not absorbed as quickly as other lifting agents:

EverLift provides margin delineation with darker color



Other Lifting Agents



EverLift

- EverLift is specially formulated with water retaining properties to prevent rapid absorption by the body

- EverLift is hypertonic, and has a slower rate of absorption due to osmotic resistance

EverLift holds a 2mm lift for 60 minutes

THE LIFT you've been waiting for



CONVENIENTLY PACKAGED AS EASY TO USE AS SPOT EX



PRE-MIXED WITH METHYLENE BLUE

Methylene blue is used to highlight the lesion perimeter and aids in:

- Increased occurrence of en bloc resection rates
- Reduced occurrence of residual lesions and incomplete resections
- Improved identification of muscle injury to help avoid missed perforations



REQUIRES NO PRE-PREP

- EverLift's convenient zero-prep, pre-filled syringe attaches directly to any endoscopic injection needle's luer-lock connection (23G recommended)
- No refrigeration required



FOR EVERYDAY USE

- Use EverLift for all polyps requiring a lift – including smaller polyps (<20mm)
- Use only what you need without throwing out what you've paid for

**First 5mL
Lifting Agent!**





For 20 years, GI Supply has been the manufacturer of the leading endoscopic tattoo, Spot Ex, which has been used by physicians and surgeons in more than 5 million cases worldwide.

GI Supply continues to broaden its submucosal injection offering with innovative and affordable devices designed to maximise procedural efficiency and patient safety.



SUPERIOR CUSHION



CONVENIENTLY PACKAGED



FOR EVERYDAY USE

Scan to learn more!



Why choose EverLift?

“This is an effective lifting agent for endoscopic resection that is offered in convenient and cost-effective packaging.”

Dr. Douglas Rex

Director of Endoscopy at
Indiana University Hospital
in Indianapolis



ORDER NO.	DESCRIPTION	QUANTITY
GIS-55	EverLift Submucosal Lifting Agent 5mL	Box of 10 Syringes

Disclaimer: EverLift Submucosal Lifting Agent 5mL is FDA approved. It is currently unavailable for sale in the EU.

For More Information Or To Locate A Distributor
Visit: www.gi-supply.com

